

— PHYSIOTHERAPY · DEEP THERMAL THERAPY

TECAR Therapy

Transfer of Energy — Capacitive & Resistive. A modern, non-invasive technique for muscular, joint and skeletal recovery.

ZE

Zahra Etminan
DPT Student · Pelvic Floor Physiotherapist




— THE ORIGIN

History

The founder of this current is Jacques Arsène d'Arsonval, a famous doctor, but also a French physicist and inventor.

In 1890, Arsonval discovered the benefits of the permeability of plasma membrane by increasing frequency to more than 100 KHz. In the 1920s, with other researchers and doctors, he developed the first capacitive and resistive electrodes.



Jacques-Arsène d'Arsonval
Physician · Physicist · Inventor

DEFINITION

TECAR

Transferring Energy CAPacitive & RESistive

100+

KHZ FREQUENCY

0

INCISIONS

What Is TECAR Therapy?

TECAR Therapy is a modern, non-invasive physiotherapy technique that uses **high-frequency electromagnetic waves** to treat muscular, joint, and skeletal injuries.

It's a form of deep thermal therapy that enhances blood circulation and stimulates the body's natural healing process by increasing internal tissue temperature.

— TWO MODES

Capacitive & Resistive

C

Capacitive Mode

TARGETS SOFT TISSUES — MUSCLES & TENDONS

Helps rapidly improve localised circulation, reduce muscle tension, and ease acute discomfort. Particularly effective for surface-level injuries and tissues that benefit from gentle, circulatory stimulation before deeper treatment is applied.

R

Resistive Mode

REACHES DEEP TISSUES — JOINTS, LIGAMENTS, BONE

CRET therapy in resistive mode generates a stronger thermal effect, stimulating cellular repair and reducing deep-seated inflammation. Commonly used for chronic joint conditions, ligament injuries, and musculoskeletal problems unresponsive to surface treatment alone.

MECHANISM

How It Works

1

Radiofrequency

Utilizes radiofrequency waves between 300 kHz and 1.2 MHz

2

Deep Penetration

These waves penetrate deep into tissues, generating internal heat

3

Healing Response

Boosts cellular metabolism, reduces inflammation, and accelerates tissue repair

CLINICAL USE

Therapeutic Applications

01 Sports injuries such as muscle and ligament strains

02 Chronic pain relief including back and neck pain

03 Arthritis and joint inflammation

04 Post-surgical rehabilitation

05 Carpal tunnel syndrome treatment

06 Tendonitis and tendon inflammation

07 Trigger point therapy

08 Reduction of swelling and lymphatic edema

09 Local emaciation

10 Infertility

**Non-invasive & painless**

No incisions, no discomfort during treatment

**No medication or surgery**

A drug-free, surgery-free recovery pathway

**Speeds up recovery**

Accelerates the body's natural healing process

**Ideal for athletes**

Trusted by professional sports performers

**Combines easily**

Works alongside other physiotherapy methods

WHY TECAR

Benefits of TECAR Therapy

A gentle yet powerful approach that brings measurable recovery without the burden of invasive procedures.

— SAFETY FIRST

Contraindications

TECAR therapy should be avoided in the following cases:



Pregnancy — particularly abdominal & pelvic regions



Pacemakers or other implanted electronic devices



Over malignant tumors or cancerous areas



Active infections, especially in the treatment area



Thrombophlebitis or deep vein thrombosis (DVT)



Severe cardiovascular disorders



Sensitivity / allergy to electromagnetic fields



Metal, cement or plastic implants & piercings at the site



Orthodontic devices (for facial-area treatment)

4—6

SESSIONS FOR ACUTE INJURIES

TREATMENT PLAN

How many sessions are usually required?

The number of sessions depends on the condition and its severity. Acute injuries may respond within **four to six sessions**, while chronic conditions typically benefit from a longer course.


PROTOCOL · CONSTIPATION

TECAR Therapy for Constipation

1 STEP

DURATION	2 minutes (deep set)
INTENSITY	40 – 60%
LOCATION	Cannon-Böhm point

Orte des Cannon-Böhm-Punkts



Cannon-Böhm point
Medi-Ko

PROTOCOL · CONSTIPATION

TECAR Therapy for Constipation

2 STEP

DURATION	Each point 2 minutes
INTENSITY	40 – 60%
LOCATION	Angles of colon




Angles of the colon

— PROTOCOL · CONSTIPATION

TECAR Therapy for Constipation

3 STEP

DURATION 15 minutes
TECHNIQUE Watch / continuous
LOCATION Colon — clockwise massage



— SPECIALISED APPLICATION

TECAR Therapy for Infertility


A targeted thermal protocol applied across the **ovary**,
fallopian tube and **uterus** regions.



PROTOCOL · INFERTILITY

How to Use

LOCATION	ASIS 2" Med./1" Inf. = ovary · 1" Med./1" Inf. = fallopian tube · 2" Med./2" Inf. = uterus
DEPTH	Maximum
INTENSITY	Generates mild, comfortable warmth
TIME	21 minutes — 7 minutes each area



ANATOMICAL REFERENCE


Female Reproductive Anatomy

The infertility protocol targets three key regions shown in the diagram — the **ovary**, **fallopian tube**, and **uterus** — applying controlled thermal energy to each area to support circulation and tissue function.

Ovary

Tube

Uterus



— QUICK REFERENCE

Treatment Parameters by Condition

Condition	Mode / Goal	TECAR Mode	Duration	Frequency	Thermal Note
Acute (Inflammatory)	Pain relief, reduce inflammation	Capacitive	10–15 min	300–500 kHz	Low (non-thermal)
Subacute	Stimulate circulation, healing	Capacitive / Resistive	15–20 min	500–800 kHz	Moderate (warm)
Chronic	Deep tissue regeneration	Resistive	20–30 min	800 kHz – 1 MHz	High (thermal)
Functional Recovery	Muscle activation, mobility	Alternating modes	15–25 min	Variable	Patient-tolerated

Parameters are indicative and should be adapted to each patient's tolerance and clinical presentation.

— TECAR THERAPY

Thank You

For your attention — questions are warmly welcomed.



Zahra Etminan

DPT Student · Pelvic Floor Physiotherapist

